



BECKER

Product range:
Endobox

INFORMATION

3-2
EN

1-2017



Magnetic endoscope holder (EBS 8-C.51)



All the components of our endoscopy product range fit together perfectly. This makes work much easier: our practical Endoboxes and manoeuvrable, individually configurable equipment trolleys create mobility and ensure that instruments are transported safely and are properly placed ready for use. Magnetic holders make sure that endoscopes are stored safely and placed in the correct positions as needed.



Seals are always at hand to mark used and unused endoscopes.



Can be machine-cleaned and disinfected, withstands temperatures up to 115°C



Endobox, grey

500 x 535 mm, 133 mm high
 W3086 70010

Lid for Endobox, grey

500 x 535 mm
 W3086 70020

Seal, package containing 500 seals
 red W9060 71050
 green W9060 72050

For other products please see www.beckerdirekt.de



131152 with Endoboxes/
lids/seals



131300 with Endoboxes/
lids/seals



Endobox equipment trolley

– suitable for pushing under other furniture –
4 pairs of holding brackets for Endoboxes/Endobox-lids

4 swivel castors 75 mm ø,
all with brakes;
Dim.: 630 x 630 mm,
height 790 mm
□ 131150

4 swivel castors 125 mm ø,
all with brakes;
Dim.: 630 x 630 mm,
height 855 mm
□ 131152

Endobox equipment trolley

5 pairs of holding brackets for Endoboxes/Endobox-lids

4 swivel castors 75 mm ø,
all with brakes;
Dim.: 630 x 630 mm,
height 955 mm
□ 131160

4 swivel castors 125 mm ø,
all with brakes;
Dim.: 630 x 630 mm,
height 1020 mm
□ 131162

Endobox transportation trolley

8 pairs of holding brackets for Endoboxes/Endobox-lids

4 swivel castors 75 mm ø,
all with brakes;
Dim.: 630 x 630 mm,
height 1560 mm
□ 131180

4 swivel castors 125 mm ø,
all with brakes;
Dim.: 630 x 630 mm,
height 1625 mm
□ 131182


Endobox transportation trolley, enclosed

pair of hinged doors, lockable, 2 push handles,
6 pairs of holding brackets for Endoboxes/Endobox-lids
4 swivel castors 125 mm ø, 2 with brakes; wall protection strips
[685 x 715 mm, height 1465 mm] □ 131300

**Endobox transportation trolley, enclosed
suitable for standard trolley washing plants** (tedesco CWA/WWA)

pair of hinged doors, lockable, 2 push handles,
6 pairs of holding brackets for Endoboxes/Endobox-lids
4 swivel castors 160 mm ø with CrNi-Steel casings,
2 with brakes; wall protection strips
[685 x 715 mm, height 1505 mm] □ 131301

*Material:
These trolleys are
made of easy-to-clean
and highly-hygienic
chromium-nickel steel,
AISI 304 (1.4301)*



Pass-through
cabinet



Cleaning-table
system with Becker
endoscope sink,
dosing device and
hygiene rail

BECKER



For decades Becker quality steel furniture and equipment have been a fixed item in medical technology in clinics and hospitals.

Design, individual choice of colours as well as an extensive range of items and sizes allow you to make optimum use of your available space. The resulting concept is customized to your needs and logical through and through.

Metallwarenfabrik
Walter H. Becker GmbH
Anzenkirchener Straße 4
84371 Triftern, Germany

Phone +49(0)8562 9605-0
Fax +49(0)8562 9605-64
info@becker-triftern.de
www.becker-triftern.de

Becker equipment:

- for the ward
- for the operating theatre
- for the central sterilisation unit
- for the endoscopy
- for all function rooms